OrthoXel

Cube House,

Model Farm Road,

Cork, T12 H966,Ireland,

+353 21 2429 500

CK-DHF1-IFU-001 Rev 11

Date: 11th July 2023

ENApex Nailing System
Attention Operating Surgeon

IMPORTANT NOTES
Intramedullary nails are intended to aid fracture healing and not to replace normal skeletal structures. Appropriate selection of implant size, length, and interlocking is essential for successful fracture treatment. As with all fracture fixation devices, healing is not guaranteed and outcomes vary by patient. Weight bearing on bones that have failed to heal or have partially or improperly healed can produce stresses in metallic implants that may lead to fatigue failure or breakage of the implant. Follow-up care and consultation with the patient are essential for safe and effective use. Metallic surgical implants are never to be reused (single use).

SYSTEM DESCRIPTION
The Apex Nailing System consists of nails of various lengths and diameters with distal and proximal holes to accept self-tapping locking screws. The Apex Femoral nail is a single product allowing implantation from an antegrade or retrograde orientation. Refer to the Apex Tibial or Femoral Nail Surgical Technique guide for a full list of sizes. The proximal end of every nail contains a preassembled insert, which is temporarily aligned and retained by means of a removable pin. The nail is designed to accept a proximal threaded end cap, which may be chosen from a range of available extension lengths. The nails, inserts, screws, and end caps are all manufactured from titanium alloy (Ti-6Al-4V EU) in conformance with ASTM and/or ISO standards.

EUROPE INDICATIONS
The Apex Nailing System is intended for temporary fixation and stabilization of the fractures of the tibia or femur. The Apex Nail is indicated for used in adult patients for treatment of:

- Open and closed tibial or femoral fractures.
- nonunion, and malunion;
- In addition the Apex Femoral Nail can treat;
- Supracondylar fractures (retrograde only);

UNITED STATES INDICATIONS
The Apex Nailing System is intended for temporary fixation and stabilization of the fractures of the tibia or femur. The Apex Nail is indicated for used in adult patients for treatment of:

- Open and closed tibial or femoral fractures, including simple, severely comminuted, spiral, large oblique, and segmental fractures;
- Fractures involving osteopenic or osteoporotic bone;
- Fractures with bone loss;
- Treatment of pseudoarthrosis, nonunion, and malunion;
- Correction osteotomy;
- Pathologic fractures and prophylactic nailing of impending pathologic fractures; In addition the Apex Femoral Nail can treat;
- Periprosthetic fractures;
- Closed supracondylar fractures (retrograde only) and;
- Intertrochanteric and subtrochanteric fractures (antegrade only).

- CONTRAINDICATIONS**
Surgeons are expected to exercise sound judgement in the selection of patients who may be suitable for intramedullary nailing. Possible contraindications for the use of this device may include:
- Active or previous infection;
 - Fixation across open epiphyseal plates;
 - Insufficient quantity or quality of bone, obliterated medullary canal, and congenital or acquired bony deformity;
 - Foreign body sensitivity including suspected or documented metal allergy or intolerance;
 - Patient conditions including limited blood supply to the injured limb, pulmonary insufficiency, hypovolemia, hypothermia, or coagulopathy;
 - Patients with mental or neurologic conditions who are unwilling or unable to follow postoperative care instructions.

- WARNINGS**
Patient Selection and Information
- The general principles of patient selection and sound surgical judgment apply to the intramedullary nailing procedure.
 - Successful fracture management requires that the patient be willing and able to adhere to the follow-up care guidelines proscribed by the surgeon. Prior to surgery, patients should be informed of the general surgical risks and possible adverse effects associated with intramedullary fixation and the expectations for follow-up examination and care.
 - Smoking increases the risk of delayed healing and nonunion. Patients should be informed of the risks associated with tobacco use and nicotine consumption in any form.
 - Nutritional insufficiency may increase the risk of delayed healing and nonunion. Patients should be counselled on the importance of a balanced and nutritious diet.
 - Patients who are senile, mentally ill, intellectually disabled, or who have uncontrolled alcoholism or substance abuse are at higher risk for non-compliance with postoperative care instructions. These patients may ignore, forget, or misunderstand activity restrictions or smoking and dietary guidelines and should be treated with judicious caution and a proactive management strategy involving all caregivers.

- Preoperative Planning**
- An adequate inventory of implant sizes, including nails, screws, and endcaps should be available at the time of surgery.
 - The Apex Nail is implanted with the Apex Nail Instrumentation Kit (as detailed in the Surgical Technique Guide), additional general surgical instruments, an intramedullary reaming system, and a surgical suite with image intensifier are required to complete the surgical procedure.
 - Test for and rule out metal allergies and sensitivities preoperatively.
 - Carefully consider the location of the fracture relative to the positions of the proximal and distal locking holes in the nail. Do not leave empty fixation holes or within the fracture line as this may increase the risk of bending or breakage of the implant.
 - When selecting an appropriate implant size for the patient, consider the fatigue life of the implant. Fatigue failure results from cyclic loading of metals over a period of time. The fatigue life (i.e. number of loading cycles before failure) is determined by the size of the nail and screws used and the loads applied by the patient according to body mass and activity level. Under equivalent loading conditions, larger nails have a longer fatigue life than smaller nails. Whenever possible, use larger nail to maximize the fatigue life of the implant system.
 - Intramedullary reaming reduces the driving force required to insert the nail and the risk of damage to the implant and the bone during nail insertion. Consider reaming 0.5-mm to 1.5-mm larger than the nail diameter chosen for the patient.

- Intraoperative Precautions**
- Surgeons are expected to have a thorough understanding of the operative procedure, instrumentation kit, and complete implant system before undertaking any procedure.
 - A purpose-built implant instrument kit has been supplied for use with this

- intramedullary nail only. The instruments are specially-designed to ensure accurate alignment of the implant components and targeted proximal drilling and screw insertion. Do not modify or supplement the instrument kit with components from any other implant system. Improper use of instrumentation may result in incorrect sizing, inaccurate drilling, operative delay, implant damage or loss of designed functionality, excessive wear, and/or premature device failure.
- A stable construct should be achieved and verified under image intensification.
 - Remove all packaging materials and instrument components prior to closing the surgical site.
 - In certain cases, a bone graft may be appropriate.

Postoperative Precautions

- Intramedullary nails are not intended to support significant weight bearing in the absence of bony union for an extended period of time. As with any fixation device, excessive cyclic loading of the metallic components causes repetitive stresses, which can lead to component fatigue fracture (i.e. breakage of screws or nail). For this reason, clear communication of appropriate early weight-bearing guidelines to the patient is critically important. Patients must be admonished that the implant system is not a replacement for normal healthy bone and that weight bearing must be limited until union can be confirmed by radiographic and clinical examination.
- Overweight and obese patients and patients who are non-compliant with regard to early weight-bearing guidelines must be advised of the risk of implant fatigue failure and be offered auxiliary support (e.g. external supports, braces, walking aids, etc.).
- Patient stature and the size and shape of the fractured bone may limit the size and strength of the implant used. Surgeons must exercise sound judgment in prescribing postoperative care and weight-bearing guidelines for all patients and especially those who require smaller-diameter nails and screws.
- Patients who are pre-disposed for healing complications or who experience a delayed union or nonunion should be cautioned against excessive weight-bearing and advised of their increased risk of implant failure due to the absence of bony support.
- Additional weight-bearing restrictions may be appropriate for patients with a fracture line occurring within 5 cm of a screw hole because this scenario may place greater stress on the nail at the location of the hole and increase the risk of implant fatigue failure.
- Patients should receive periodic clinical and radiographic evaluations for at least the first six months post-fixation. Monitor patients for evidence of implant migration, loosening, bending, or breakage. If any of these conditions arise, especially in cases of delayed or nonunion, implant damage may worsen, so consider additional activity restrictions or early revision.
- Diffcult fractures may take 8-12 months to heal. Patients who experience a delayed union are at increased risk for implant fatigue failure, especially when smaller diameter nails are used. For high-risk patients who require a small-diameter nail, consider one or more of the following adjunctive treatments: fracture-site grafting at the time of surgery, or prophylactic exchange nailing for a larger-diameter nail at 3-4 months if there is no radiographic evidence of union.
- After complete union has been achieved, the presence of the implant may cause stress shielding, which may increase the risk of re-fracture in an active patient. The surgeon should weigh the risks versus benefits when deciding whether to remove the implant after healing is complete, especially in elderly or debilitated patients.
- If implant removal is indicated, provide the patient with appropriate postoperative management guidelines to reduce the risk of re-fracture after implant removal. Counsel the patient that this re-fracture is most likely soon after explantation, when voids left in the bone after implant removal have not filled in completely.
- After removal, label all Implant components as a potential biohazards. Reuse of devices labelled for single-use may result in product contamination and patient infection.

Other Mechanical and Metallurgical Considerations

- Surgeons are expected to understand the medical, surgical, mechanical, and metallurgical properties, features, and limitations of the implant system.
- In addition to the risk of fatigue failure under excessive loading conditions, implants may fail if improperly handled pre- or intra-operatively. Do not bend, notch, scratch, cut, or otherwise modify nails or screws. Improper handling or modification may increase the risk of implant breakage.
- Do not apply excessive torque when seating bone screws; this may cause intraoperative screw breakage.
- Implanted metals and metal alloys are subject to corrosion due to constantly-changing in vivo conditions that may expose the implant to salts, acids, and alkalis.
- The surfaces of this implant system have been treated by anodisation to reduce their corrosion potential. Scratching, filing, notching, abrading, or otherwise removing or disturbing the surface coating may increase the risk of corrosion and must be avoided.
- Dissimilar metals and metal alloys in contact with each other are at increased risk of corrosion. All components of this nailing system are made from the same metal.
- Paints or finishes should be combined with elements from another nailing system.
- All components of the implant system (nail, insert, screws, and endcap) are single-use only. Never reuse any component of the implant system, even if it was implanted only momentarily and has been subsequently sterilized. A previously-used implant may appear to be undamaged, but stresses generated during the first implantation may have caused damage that could reduce the service life of the implant.

POSSIBLE ADVERSE EFFECTS

- Nonunion or delayed union.
- Malunion with loss of anatomic alignment.
- Limb shortening due to compression of the fracture or bone resorption.
- Loosening or migration of the implant, including loss of fixation due to osteoporosis.
- Bending or breakage of the implant.
- Scar formation and/or soft tissue injury.
- Nerve damage due to surgical or pre-existing trauma.
- Pain, discomfort, or abnormal sensation due to the presence of the implant.
- Metal sensitivity or allergic reaction to a foreign body.
- Increase in bone density due to stress shielding.
- Necrosis of bone.
- Postoperative bone fracture and/or pain.
- Inadequate healing.
- Early or late infections, both deep and superficial.

MRI INFORMATION

Non-clinical testing has demonstrated that the Apex Nailing System is MR Conditional: A patient with this device can be safely scanned in an MR system meeting the following conditions:

Cylindrical-bore,
Horizontal magnetic field of 1.5 or 3 T,
Maximum spatial field gradient: (1.5 T: 73.64 T/m, 3 T: 36.82 T/m),
Radio-frequency (RF) field exposure, RF excitation:
Circularly Polarized (CP),
RF transmit coil type: whole-body transmit coil,
RF receive coil type: whole-body receive coil,
Maximum permitted whole-body averaged specific absorption rate (SAR):
Normal Operating Mode,
2 W/kg,
Scan duration: Maximum 1h of continuous scanning.
Metal implants can result in image artifacts that could cause misinterpretation of results. When deciding whether patients should undergo MRI, physicians must consider both the advantages of imaging and the possibility of image distortion due to the implants.

STERILITY

OrthoXel intramedullary Nail implants are provided sterile and clearly marked "Sterilized using irradiation" on the packaging. These implants have been sterilized using a minimum dosage of 2.5 megaRad (25 kGy) of gamma radiation. Where specified, do not use implants after expiration date.

PACKAGING

Store packaging in a cool dry location away from direct sunlight. If sterile packaging is damaged do not use implant.

CAUTION:

Federal law (USA) restricts this device to sale by or on the order of a physician.

INFORMATION

Comments regarding this device can be directed to Attn: Regulatory Dept., OrthoXel, Cube House, Model Farm Road, Cork, T12 H966, Ireland, Tel: + 353 21 2429 500
All trademarks herein are the property of OrthoXel or its subsidiaries unless otherwise indicated. CE Mark on the package insert (IFU) is not valid unless there is a CE Mark on the product (description) label.

REPORTING SERIOUS INCIDENTS (EU)

An y serious incident that has occurred in relation to this system should be reported to OrthoXel and the competent authority of the Member State in which the reporting party is established.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (EU)

A summary of safety and clinical performance can be found at www.orthoXel.com/resources/#sscp

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