EN **OrthoXel: Locking Screws and Bolts**
Attention Operating Surgeon

IMPORTANT NOTES

Intramedullary nails and the associated locking screws and bolts are intended to aid fracture healing and not to replace normal skeletal structures. Appropriate selection of implant size, length, and interlocking is essential for successful fracture treatment. As with all fracture fixation devices, healing is not guaranteed and outcomes vary by patient. Weight bearing on bones that have failed to heal or have partially or improperly healed can produce stresses in metallic implants that may lead to fatigue failure or breakage of the implant. Follow-up care and consultation with the patient are essential for safe and effective use. Metallic surgical implants are never to be reused (single use).

SYSTEM DESCRIPTION

OrthoXel: Locking Screws and Bolts have been designed to fit with the various intramedullary nail sizes. For full details, refer to the appropriate surgical technique guide for the chosen intramedullary nail. Both the screws and bolts are manufactured from titanium alloy (Ti-6Al-4V ELI) in conformance with ASTM and/or ISO standards.

INDICATIONS

Locking Screws and Bolts together with the appropriate intramedullary nail are intended for temporary fixation and stabilization of tibial or femoral fractures. Refer to the intramedullary nail IFU for the specific indications associated with the relevant nailing system.

CONTRAINDICATIONS

Refer to the appropriate intramedullary nail IFU for the specific contraindications associated with the relevant nailing system.

WARNINGS

Patient Selection and Information

- The general principles of patient selection and sound surgical judgment apply to the intramedullary nailing procedure.
- Successful fracture management requires that the patient be willing and able to adhere to the follow-up care guidelines proscribed by the surgeon. Prior to surgery, patients should be informed of the general surgical risks and possible adverse effects associated with intramedullary fixation and the expectations for follow-up examination and care.
- Smoking increases the risk of delayed healing and nonunion. Patients should be informed of the risks associated with tobacco use and nicotine consumption in any form.
- Nutritional insufficiency may increase the risk of delayed healing and nonunion. Patients should be counselled on the importance of a balanced and nutritious diet.
- Patients who are senile, mentally ill, intellectually disabled, or who have uncontrolled alcoholism or substance abuse are at higher risk for non-compliance with postoperative care instructions. These patients may ignore, forget, or misunderstand activity restrictions or smoking and dietary guidelines and should be treated with judicious caution and a proactive management strategy involving all caregivers.

Preoperative Planning

- An adequate inventory of implant sizes, including nails, screws, and endcaps should be available at the time of surgery.
- Additional general surgical instruments and a surgical suite with image intensifier are required to complete the surgical procedure.
- Test for and rule out metal allergies and sensitivities preoperatively.
- Carefully consider the location of the fracture relative to the positions of the proximal and distal locking holes in the nail. Do not leave empty fixation holes at or within the fracture line as this may increase the risk of bending or breakage of the implant.
- When selecting an appropriate implant size for the patient, consider the fatigue life of the implant. Fatigue failure results from cyclic loading of metals over a period of time. The fatigue life (i.e. number of loading cycles before failure) is determined by the size of the nail and screws used and the loads applied by the patient according to body mass and activity level. Under equivalent loading conditions, larger nails, with larger screws have a longer fatigue life.
- Intramedullary reaming reduces the driving force required to insert the nail and the risk of damage to the implant and the bone during nail insertion. Consider reaming 0.5-mm to 1.5-mm larger than the nail diameter chosen for the patient.

Intraoperative Precautions

- Surgeons are expected to have a thorough understanding of the operative procedure, instrumentation kit, and complete implant system before undertaking any procedure. Refer to the Surgical Technique for a complete system description and procedural notes.
- A purpose-built implant instrument kit has been supplied for use with the OrthoXel: Locking Screw and Bolts only. The instruments are specially designed to ensure accurate alignment of the implant components and targeted proximal drilling and screw insertion. Do not modify or supplement the instrument kits with components from any other implant system. Improper use of instrumentation may result in incorrect sizing, inaccurate drilling, operative delay, implant damage or loss of designed functionality, excessive wear, and/or premature device failure.
- A stable construct should be achieved and verified under image intensification.
- Remove all packaging materials and instrument components prior to closing the surgical site.
- In certain cases, a bone graft may be appropriate.

Postoperative Precautions

- Intramedullary nails and the associated locking screws and bolts are not intended to support significant weight bearing in the absence of bony union for an extended period of time. As with any fixation device, excessive cyclic

loading of the metallic components causes repetitive stresses, which can lead to component fatigue fracture (i.e. breakage of screws or nail). For this reason, clear communication of appropriate early weight-bearing guidelines to the patient is critically important. Patients must be admonished that the implant system is not a replacement for normal healthy bone and that weight bearing must be limited until union can be confirmed by radiographic and clinical examination.

- Overweight and obese patients and patients who are non-compliant with regard to early weight-bearing guidelines must be advised of the risk of implant fatigue failure and be offered auxiliary support (e.g. external supports, braces, walking aids, etc.).
- Patient stature and the size and shape of the fractured bone may limit the size and strength of the implant used. Surgeons must exercise sound judgement in proscribing postoperative care and weight-bearing guidelines for all patients and especially those who require smaller-diameter nails and screws.
- Patients who are pre-disposed for healing complications or who experience a delayed union or nonunion should be cautioned against excessive weight-bearing and advised of their increased risk of implant failure due to the absence of bony support.
- Additional weight-bearing restrictions may be appropriate for patients with a fracture line occurring within 5 cm of a screw hole because this scenario may place greater stress on the nail at the location of the hole and increase the risk of implant fatigue failure.
- Patients should receive periodic clinical and radiographic evaluations for at least the first six months post-fixation. Monitor patients for evidence of implant migration, loosening, bending, or breakage. If any of these conditions arise, especially in cases of delayed or nonunion, implant damage may worsen, so consider additional activity restrictions or early revision.

- Difficult fractures may take 8-12 months to heal. Patients who experience a delayed union are at increased risk for implant fatigue failure, especially when smaller nails are used. For high-risk patients who require a small-diameter nail, consider one or more of the following adjunctive treatments: fracture-site grafting at the time of surgery, nail dynamization at 2-3 months, or prophylactic exchange nailing for a larger-diameter nail at 3-4 months if there is no radiographic evidence of union.
- After complete union has been achieved, the presence of the implant may cause stress shielding, which may increase the risk of re-fracture in an active patient. The surgeon should weigh the risks versus benefits when deciding whether to remove the implant after healing is complete, especially in elderly or debilitated patients.
- If implant removal is indicated, provide the patient with appropriate postoperative management guidelines to reduce the risk of re-fracture after implant removal. Counsel the patient that this re-fracture is most likely soon after explantation, when voids left in the bone after implant removal have not filled in completely.

- After removal, label all implant components as a potential biohazards. Reuse of devices labelled for single use may result in product contamination and patient infection.

Other Mechanical and Metallurgical Considerations

- Surgeons are expected to understand the medical, surgical, mechanical, and metallurgical properties, features, and limitations of the implant system.
- In addition to the risk of fatigue failure under excessive loading conditions, implants may fail if improperly handled pre- or intraoperatively. Do not bend, notch, scratch, cut, or otherwise modify nails or screws. Improper handling or modification may increase the risk of implant breakage.
- Do not apply excessive torque when seating bone screws; this may cause intramediated screw breakage.
- Implanted metals and metal alloys are subject to corrosion due to constantly-changing in vivo conditions that may expose the implant to salts, acids, and alkalis.
- The surfaces of this implant system have been treated by anodisation to reduce their corrosion potential. Scratching, filing, notching, abrading, or otherwise removing or disturbing the surface coating may increase the risk of corrosion and must be avoided.
- Dissimilar metals and metal alloys in contact with each other are at increased risk of corrosion. All components of this nailing system are made from the same metal alloy and should never be combined with elements from another nailing system.
- All components of the implant system are single-use only. Never reuse any component of the implant system, even if it was implanted only momentarily and has been subsequently sterilized. A previously used implant may appear to be undamaged, but stresses generated during the first implantation and explantation may have caused damage that could reduce the service life of the implant, change the functional performance of the device, and/or cause harm if used again.

POSSIBLE ADVERSE EFFECTS

- Nonunion or delayed union.
- Malunion with loss of anatomic alignment.
- Limb shortening due to compression of the fracture or bone resorption.
- Loosening or migration of the implant, including loss of fixation due to osteoporosis.
- Bending or breakage of the implant.
- Scar formation and/or soft tissue injury
- Nerve damage due to surgical or pre-existing trauma.
- Pain, discomfort, or abnormal sensation due to the presence of the implant.
- Metal sensitivity or allergic reaction to a foreign body.
- Decrease in bone density due to stress shielding.
- Necrosis of bone.
- Postoperative bone fracture and/or pain.
- Inadequate healing.
- Early or late infections, both deep and superficial.
- Where applicable, Screw cutout of the femoral head (usually associated with osteoporotic bone).

MRI SAFETY INFORMATION:

Refer to the appropriate intramedullary nail IFU for the specific MRI Safety Information associated with the relevant nailing system.

STERILITY

OrthoXel implants are provided sterile and clearly marked “Sterilized using irradiation” on the packaging. These implants have been sterilized using a minimum dosage of 2.5 megaRad (25 kGy) of gamma radiation. Where specified, do not use implants after expiration date.

PACKAGING

Store packaging in a cool dry location away from direct sunlight. If sterile packaging is damaged, do not use implant.

CAUTION:

Federal law (USA) restricts this device to sale by or on the order of a physician.

INFORMATION

Comments regarding this device can be directed to Attn: Regulatory Dept., OrthoXel, Cube House, Model Farm Road, Cork, T12 H966, Ireland. Tel: + 353 21 2429 500 All trademarks herein are the property of OrthoXel or its subsidiaries unless otherwise indicated.

CE Mark on the package insert (IFU) is not valid unless there is a CE Mark on the product (description) label.

REPORTING SERIOUS INCIDENTS (EU)

Any serious incident that has occurred in relation to this system should be reported to OrthoXel and the competent authority of the Member State in which the reporting party is established.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (EU)

A summary of safety and clinical performance can be found at www.orthoxel.com/ resources/#sscsp

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Symbol Key
Legenda
Légende
Leyenda de símbolos utilizados
Zeichenerklärung

	Manufactured in accordance with the requirements of the Medical Device Regulation (MDR) 2017/745.		Read the instructions for use carefully before use.
	Do not smoke during the operation.		Do not use open flames or heat sources near the device.
Qty	Quantity		Do not use the device in wet conditions.
REF	Reference number		Do not use the device if you are under the influence of alcohol or drugs.
Serial R	Serial number		Do not use sharp objects to clean the device.
LOT	Lot number		Do not use sharp objects to clean the device.
	Do not smoke during the operation.		Do not use sharp objects to clean the device.
	Do not smoke during the operation.		Do not use sharp objects to clean the device.
	Do not smoke during the operation.		Do not use sharp objects to clean the device.
	Do not smoke during the operation.		Do not use sharp objects to clean the device.
	Do not smoke during the operation.		Do not use sharp objects to clean the device.